

PROHS Equipamento Hospitalar
e Serviços Associados, S.A.

Rua do Castanhal 316
Zona Industrial Maia I, sector II
4475-122 Maia, Portugal

| T | +351 229 059 170
| F | +351 229 015 900
| E | prohs@prohs.pt
| W | www.prohs.pt

Closed Car for Sterilized Items

PROHS[®]



Closed Car for Sterilized Items



Model N-6

Intended for transportation and distribution of containers and /or baskets with sterilized material or return the same material to the CSSD.

Technical Characteristics

- Built in stainless steel AISI 304L.
- Interior rails for wired baskets.
- The hermetic seal achieved by placing a rubber gasket on the door.
- Breaks on both wheels.
- Removable shelves.
- Capacity for 6 (N-6 model) or 9 (N-9 model) sterilization units.
- The N-9 model is equipped with two doors allowing an opening to 270 °.



Model N-9

Brand	Model	Dimensions (mm) W. x D. x H.
PROHS	N-6	660 x 710 x 1220
	N-9	1200 x 695 x 1450

Can be produced in different dimensions.

Due to product improvement and development, PROHS reserves the right to make changes without notice.

