


# Preparation and Packing Table

PROHS Equipamento Hospitalar  
e Serviços Associados, S.A.

Rua do Castanhal, 316  
Zona Industrial da Maia I, Sector II  
4475-122 Maia, Portugal

| T | +351 229 059 170  
| F | +351 229 015 900  
| E | [prohs@prohs.pt](mailto:prohs@prohs.pt)  
| W | [www.prohs.pt](http://www.prohs.pt)



# Preparation and Packing Table

The Preparation and Packing Tables are specially developed for the sterilization rooms.



Preparation and Packing Table - double post  
2000 x 850 x 900mm

## Technical Characteristics

- Made in stainless steel AISI 304L with flat tops and soundproofed, the support consists of a welded tubular Argon structures.
- Rounded and polished finish for easy cleaning and promoting hygiene products.
- Structure of square tubing in stainless steel AISI 304, with adjustable feet height.
- The top shelf, tilted for easy access by the operator. Shelf less than the length of the table.
- Accommodate various accessories: paper holder, holder baskets, drawers, etc.
- Possibility of having two opposite positions.

**Can be produced in different dimensions.**



Preparation and Packing Table - simple  
2000 x 1000 x 900mm



Preparation and Packing Table - simple post  
1500 x 700 x 900mm

Due to product improvement and development, PROHS reserves the right to make changes without notice.

