

## PROHS SEALERS GENERAL CATALOGUE

PROHS Equipamento Hospitalar  
e Serviços Associados, S.A.

Rua do Castanhal, 316  
Zona Industrial Maia I, sector II  
4475-122 Maia — Portugal

| E | [prohs@prohs.pt](mailto:prohs@prohs.pt)  
| W | [www.prohs.pt](http://www.prohs.pt)



# PROHS SEALERS

	TS 45 DN	TS 46N	e-FLEX	TS 47N	e-PRINT
<b>TECHNICAL CHARACTERISTICS</b>					
Type: Rotary Sealer	✓	✓	✓	✓	✓
Printing (texts and symbols): stamping date, expiry date, users, batches, stamping counter, temperature, pressure, free text, ...	-	-	-	✓	✓
<b>STANDARDS</b>					
CE / CEM	✓	✓	✓	✓	✓
EN 868-5	✓	✓	✓	✓	✓
DIN 58953-7	✓	✓	✓	✓	✓
EN ISO 11607-2	✓	✓	✓	✓	✓
Temperature calibration certificate (COFRAC)	option	option	option	option	option
<b>Screen</b>					
Screen	1 line LCD	1 line LCD	7" Tablet	1 line LCD	7" Tablet
Intuitive multifunctional menu (graphs, counters,informative messages, multilingual, alphanumeric keyboard ...)	-	-	✓	-	✓
Digital maintenance notebook	-	-	✓	-	✓
Temperature range	20-230 °C	20-230 °C	20-230 °C	20-230 °C	20-230 °C
Pre-programmed temperatures	-	-	✓	-	✓
Temperature regulation tolerance	2%	2%	2%	2%	2%
<b>SAFETY</b>					
Automatic switch-on at set point	✓	✓	✓	✓	✓
Automatic stop if overheating (250° C)	✓	✓	✓	✓	✓
Control & Validation of critical parameters in accordance with ISO 11607-1 (Temperature / Force / Speed)	T° / F	T° / F	T° / F / S	✓	✓
Automatic stop if parameters exceed limits	✓	✓	✓	✓	✓
<b>OTHER INFORMATION</b>					
Automatic start and stop by photoelectric cell	✓	✓	✓	✓	✓
Adjustable bag guide	fixed 30mm	0 to 30mm	0 to 30mm	0 to 30mm	0 to 30mm
Sealing width	15mm	15mm	15mm	15mm	15mm
Sealing length	Unlimited	Unlimited	Unlimited	Unlimited	Unlimited
Types of bags: Steam/ OE / alu / Tyvek® / Ultra®	✓	✓	✓	✓	✓
Electrical / frequency connection	230V - 50/60 Hz	230V - 50/60 Hz	230V - 50/60 Hz	230V - 50/60 Hz	230V - 50/60 Hz
Power	400 W max.	400 W max.	400 W max.	400 W max.	400 W max.
Cover	Stainless steel	Stainless steel	Stainless steel	Stainless steel	Stainless steel
Dimensions - Length x Width x Height (mm)	380 x 240 x 150	439 x 244 x 162	439 x 244 x 162	555 x 265 x 178	455 x 265 x 178
Weight (kg)	8	10,5	11	13,5	16
<b>ACCESSORIES &amp; OPTIONS</b>					
Front tray in Stainless steel (620 x 210 x 50mm)	option	option	option	option	option
Roller conveyor (800 x 295 x 50mm)	option	option	option	option	option
OQ EDIT - integrated self-test printer for the editing of monitoring reports ISO 11607-2 (Temperature, Force, Speed, Maintenance)	-	-	option	-	option
Traceability and storage of all parameters via USB memory or LAN (Ethernet)	-	-	option	-	option
Barcode and scanner software for programming of sealers (SCAN)	-	-	option	✓	✓
USB Connection	-	-	option	-	-



TS 45DN with Front tray



e-FLEX with Roller conveyor



e-FLEX with USB

